

Application Brief

Blood Bag Labeling

Industry:
Healthcare

Products:
Duratran II Blood Bag Label
TMX3000 series ribbon
TMX2000 series ribbon
PM43
PC23/PC43

Region:
North America

Typical Applications:
Blood bag labeling

Customer Benefits:

- Reliable blood bag tracking with durable labels
- Assured compliance with government and industry standards
- Application-tailored adhesive and face sheet to maintain bond with flexible blood bags during processing and handling
- Durable, high contrast print for consistent scanning performance
- Reduced total cost of ownership with free printhead replacement program when Intermec media and printers are used together



Industry Need

With nearly 5 million Americans requiring blood transfusions annually, there is a tremendous need for blood products to supply these important therapies. Blood banks and collection centers require robust labels to reliably track blood from the donor to the final recipient; these labels need to survive multiple processing, testing, and storage steps through challenging environmental conditions. Intermec's blood bag labeling product set enables durable, positive tracking of these critical components.

Bar coded blood bag labels enable the high degree of accuracy required for successful transfusion therapy. Clearly identified blood typing prevents serious complications that could result from donor type incompatibility. Different blood components and storage conditions yield varying shelf life; clear expiration date identification guarantees patients receive only safe, effective products. Formats configured in compliance with ISBT-128 standards provide interoperability with multiple collection centers, processing locations, and hospital systems.

In order for essential data to follow the contents of the blood bag, labels must endure a series of challenging conditions. Labels must maintain a strong bond to the flexible blood bag during centrifugation while whole blood is separated into components. Labels must remain positively attached during standard or cryogenic freezing (plasma and some red cells), refrigeration (red cells), or continuous shaking at room temperature (platelets). The bond must then be maintained during warming, which often includes immersion in a water bath.

Beyond the physical performance needs, blood bag labels must also meet regulatory requirements. The United States Food and Drug Administration mandates compliance with 21 CFR175.105 to increase safety should components of the adhesive migrate through the plastic into the blood bag.

Reference Links

Standard for the Uniform Labeling of Cellular Therapy Products Using ISBT 128

http://icbba.org/usconsensusstandard_cellulartherapy.pdf

FDA Guideline for the Uniform Labeling of Blood and Blood Components

<http://www.fda.gov/downloads/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Blood/UCM080974.pdf>

Intermec Solution

Intermec's blood bag label products are designed for the needs of this critical tracking application. An aggressive adhesive securely bonds the label to the filled bag; the pliable face sheet flexes with it during centrifugation and multiple handling steps. With a wide service temperature range, Intermec labels remain positively attached through the broad range of blood component storage temperatures. Qualification testing through multiple freeze/thaw cycles ensures labels readily endure the harsh treatments commonly encountered.

Intermec blood bag media complies with the industry standard 21 CFR 175.105, increasing assurance of safety.

Images printed on Intermec blood bag labels with the durable TMX3000 or TMX2000 series ribbons, readily endure extensive handling, further adding to the robustness of the product set. Crisp, high contrast print ensures consistent scanning and imaging performance, eliminating error-prone manual data entry.

Labels can be provided in full face, single or double quadrant with or without DIN (Donor ID Number) cutouts. Piggyback material is also available for multi-step bag tracking or donor set labels.

Products

Duratran II Permanent Blood Bag Label

Ribbon Recommendation

TMX3000 series resin

TMX2000 series mid-range

Printer Recommendation

PM43 Mid Range Printer [Profile](#) (High volume output)

PC43t Desktop Printer [Profile](#) (Lower volume desktop use)

- **Intuitive** – *Super Easy to Install & Use*
 - 10-language LCD or Icon graphics user interface maximizes uptime
 - One hand media loading for trouble-free use
 - Fastest throughput in its class at 8 ips
 - Quick configuration with USB stick device
- **Flexible** – *Fits your business needs today & tomorrow*
 - User installable connectivity upgrades and accessories for scalability
 - Leading command languages for seamless deployment
 - Modular IT-friendly design speeds up deployment and support
- **Smart** – *Keeps you informed & productive*
 - Remote status updates with Intermec tools minimize downtime
 - Programmable with ability to directly add keyboard or scale improves efficiency

The Intermec Advantage

Intermec printers and media are optimized to deliver superior performance when used together. Our rigorous testing and co-engineering ensures maximum print head service life, consistently high print quality, and proven label and tag performance in



demanding real-world environments. Intermec printers and media are complemented by our complete line of integrated solutions and services, including rugged mobile computers, RFID systems, and advanced scanning technology. Intermec products and services are used by customers worldwide to improve the productivity, quality and responsiveness of their business operations.

Printhead Replacement Program

Intermec will provide no cost replacement printheads, as failures occur due to normal wear and tear, to those customers who use Intermec Media Products exclusively with their Intermec direct thermal and thermal transfer printers. For full details, click [here](#).

For more information

Contact Intermec for more information on robust blood bag identification.

North America

Corporate Headquarters

6001 36th Avenue West
Everett, Washington 98203
Phone: (425) 348-2600
Fax: (425) 355-9551

North Latin America

Headquarters Office

Mexico
Phone: +52 55 52-41-48-00
Fax: +52 55 52-11-81-21

South Latin America

Headquarters Office

Brazil
Phone: +55 11 3711-6776
Fax: +55 11 5502-6780

Europe/Middle East & Africa

Headquarters Office

Reading, United Kingdom
Phone: +44 118 923 0800
Fax: +44 118 923 0801

Asia Pacific

Headquarters Office

Singapore
Phone: +65 6303 2100
Fax: +65 6303 2199

Internet

www.intermec.com
Worldwide Locations:
www.intermec.com/locations

Sales

Toll Free NA: (800) 934-3163
Toll in NA: (425) 348-2726
Freephone ROW:
00 800 4488 8844
Toll ROW: +44 134 435 0296

OEM Sales

Phone: (425) 348-2762

Media Sales

Phone: (513) 874-5882

Customer Service and Support

Toll Free NA: (800) 755-5505
Toll in NA: (425) 356-1799

Copyright © 2012 Intermec Technologies Corporation. All rights reserved. Intermec is a registered trademark of Intermec Technologies Corporation. All other trademarks are the property of their respective owners. Printed in the U.S.A. 612018 04/12

In a continuing effort to improve our products, Intermec Technologies Corporation reserves the right to change specifications and features without prior notice.

